Test 3


<u>Description</u>: Estimation of the evaporation in the plate when dispensing ddH20 in 24h.

<u>**Purpose**</u>: evaluation of the homogenity of the temperature inside the instrument, find an efficient protocol of dispensation

Methods A flat-bottom non sterile plate is used. All 96 wells are filled with 200 μl LB+Kan

Protocol:

- The plate is filled as described in Methods
- The instrument temperature was set at 37°C
- Experiment duration: 24h
- Dispensations: 20µl ddH20 every 2h
- Shaking 15s linear 3mm after every dispensation.




Notes:

After the 4th dispensation (80µl ddH2O dspensed in total) the plate was took out of the reader and the final volumes in 4 wells was misured:

C7: 254 µl E12: 255 µl F10: 255 µl D1: 275 µl

A visual analisy of the plate volume's profile was done, and a disomogenity was noted.



A new dispensing protocol was chosen:

• 10 µl every 6 h

The first dispensation was done at 8.30pm, so the other will be at 2.30am and 8.30am.

The final total volume dispensed will be 80 μ l + 30 μ l : 110 μ l.