


culture dilutions TEST 07/07/09

Description: Dispensation Test

Purpose: we want to estimate how varies absorbance and fluorescence in red and green in cultures (J100, A1, RBS-30) with periodical dispensation of ddH2O

Methods: A flat-bottom non sterile plate is used. 96 wells are filled with:

- 200 ul LB+Amp,
- 200 ul bacterial cultures incubated overnight 37°C 220 rpm and then diluted 1:1000
- 150 ul bacterial cultures incubated overnight 37°C 220 rpm and then diluted 1:1000
- 100 ul bacterial cultures incubated overnight 37°C 220 rpm and then diluted 1:1000
- 50 ul bacterial cultures incubated overnight 37°C 220 rpm and then diluted 1:1000
- 200ul, 150ul, 100ul and 50ul ddH2O.



Protocol:

- The plate is filled as described in Methods
- The instrument temperature is set at 37°C
- Dispensation of 200ul, 150ul, 100ul and 50ul ddH2O in empty wells (first line, A1:a12)
- First static measures:
 - Measure of ASB, GFP, RFP in wells filled with 200ul, 150ul, 100ul, 50ul.
 - Dispensation of 50ul ddH2O in line D,E,F
 - Shaking 15 s 6mm linear

- 15s waiting
 - Measure of ASB, GFP, RFP
- Second static measures:
 - Waiting 5 minutes, shaking 15s linear 3mm, 10sec waiting
 - Measure of ASB, GFP, RFP
 - Dispensation of 50ul ddH2O in line E,F
 - Shaking 15 s 6mm linear
 - 10s waiting
- Measure of ASB, GFP, RFP
- Third static measures:
 - 5minutes waiting, 15 s shaking 3mm linear, 15 s waiting
 - Measure of ASB, GFP, RFP
 - Dispensation of 50ul ddH2O in line E
 - Shaking 15 s 3mm linear
 - 10s waiting
 - Measure of ASB, GFP, RFP